

PEN REFRACTOMETER



Cat.No. 3730

NEW Continuous

Measurement Feature

expands traditional

measuring methods



Benefits:

Fast measurements! Non-homogenous samples—no problem!

Stir the sample while taking continuous measurements until the readings become stable.

Easy Continuous Measurements!


Once the START key is pressed, measurements are taken every second.



Caution: For sampling of liquids of similar concentration and of minimum viscosity.









PEN REFRACTOMETER PFN-PK Patent Pending


Names and Functions of Components



The PEN type refractometer is extremely easy to use.










Features

Extremely water resistant (IP65) **External Light Interference (ELI) Automatic Temperature Compensation (ATC)**

Calibration with water Measurement in 2 seconds Light & Compact, 70g



External Light Interference (ELI)

When intense light passes through the prism of a digital refractometer, the light waves may interfere with the sensor, which can lead to inaccurate measurements. To ensure accurate measurement results, the PEN-PRO is programmed with the ELI function, which displays the [nnnn] warning message when intense direct light is detected. When the ELI warning is displayed, shading the prism head with your hand and pressing the START key again will ensure accurate measurements.

Specifications

Measurement range	Brix 0.0 to 85.0%
Resolution	Brix 0.1%
Measurement accuracy	Brix ±0.2%
Measurement temperature	10 to 100℃ (Automatic Temperature Compensation)
Sample temperature	10 to 100℃
Ambient temperature	10 to 40℃
Power supply	1×Size AAA alkaline battery
International Protection class	IP65 Water Resistant (Prism head:IP67)
Dimensions and weight	16(W)×3.8(D)×1.8(H)cm, 70g(main unit only)

Brix(%)

All refractometers are designed to measure the Refractive Index of a solution. The Brix scale is based on a sucrose and water solution. Since most samples contain substances other than sugar (salts, minerals, proteins), the Brix percentage represents the total concentration of all soluble solids. For certain samples, such as cutting oil or industrial fluids, a conversion chart from Brix percentage to the total concentration of the sample may be necessary. For further information, contact an ATAGO Customer Service Representative.

All ATAGO refractometers are designed and manufactured in Japan.





Headquarters: 32-10, Honcho, Itabashi-ku, Tokyo 173-0001 Japan TEL: 81-3-3964-6156 FAX: 81-3-3964-6137 overseas@atago.net http://www.atago.net/



ATAGO U.S.A., Inc.

12011 NE First Street, Bldg. C, Suite 110, Bellevue, WA 98005 U.S.A. TEL: 1-425-637-2107 FAX: 1-425-637-2110 customerservice@atago-usa.com http://www.atago.net/USA/



A-501, Mangal Aarambh Business Centre, Near Kora Kendra, Off. S.V. Road, Borivali (West), Mumbai-400 092 India TEL: 91-22-2833-8038 / 8076 FAX: 91-22-2899-8164 customerservice@atago-india.com

