

Easy to Use

High Performance

Great Low Price

# Honey Refractometer % Water Content!!

Cat. No.4422

## Digital Hand-held "Pocket" Honey Refractometer

# PAL-22S

The "Pocket" Honey Refractometer, PAL-22S can measure the % water content of honey easily and displays readings digitally. Comfortable, easy to use, light & compact. Measurements can be taken anywhere.

### Features

- \* Calibration with water only
- \* Measurement in 3 seconds
- \* Highly water resistant (IP65)
- \* Light & Compact, 100g
- \* Equipped with a strap holder option
- \* Ergonomically designed for one-handed operation
- \* Convenient storage case
- \* Automatic Temperature Compensation (ATC)



Actual size

### Measuring Method



- ① Place the sample onto the prism surface.
- ② Press the START key.



- ③ The measurement value will be displayed in 3 seconds.

### Specifications

Measurement range	Honey Moisture 12.0~30.0%
Resolution	Honey Moisture 0.1%
Measurement accuracy	Honey Moisture ±0.2%
Measurement temperature	10~40°C (ATC)
Sample volume	0.3ml or more
Measuring time	3 seconds
Power supply	AAA alkaline battery x 2
Battery life	Approx. 11,000 measurements (Using a new alkaline battery)
International Protection Class	IP65
Dimensions and weight	55(W) x 31(D) x 109(H)mm, 100g

CE \*Specifications and appearance are subject to change without notice.

**ATAGO CO., LTD.**  
 Main Office: 32-10, Honcho, Itabashi-ku, Tokyo 173-0001, Japan  
 TEL: 81-3-3964-6156 FAX: 81-3-3964-6137  
 export@atago.net http://www.atago.net/

**ATAGO U.S.A., Inc.**  
 13005 NE 126th Place, Kirkland, WA 98034 U.S.A.  
 TEL: 1-425-821-9050 FAX: 1-425-825-5525  
 customerservice@atago-usa.com http://www.atago.net/



Printed in Japan